# ID N° 20326 - Contact Angle Measurement - M4<sup>s+</sup>

## Isooktan Standard-Wilhelmy (4mm Immersion)

• **Wilhelmy-Platte (DIN) / Isooktan 99.5+%** : (20,0°C, 4,3', 4,0mm, 0,0286mm/s - static;  $\theta$ -superwetting)

$\theta_{C,s} = 0 \pm 0^\circ$  CAH  $0^\circ$ ,  $^a\bar{E}_s = 18,887 \pm 0,021$  mN/m  $^aH_s = -0,039$  mN/m,  $\bar{B}_{\%,s} = 100,7\%$

## Report

### 1. ===== Introduction =====

**Principle:** The M4 contact angle measurement using the *modified Wilhelmy method* is based on force measurements on a cuboid or cylindrical solid body brought into contact with a liquid surface. The contact line between the solid, the liquid surface, and the gas phase is called the three-phase boundary line (triple line), and the angle at which the liquid surface meets the three-phase boundary line is the contact angle  $\theta$  (theta). The relative position of the three-phase boundary line-above, on, or below the surrounding liquid level-reflects the surface energy characteristics and is indicated by the resulting curvature of the meniscus. The weight of this meniscus curvature is determined gravimetrically, taking buoyancy forces into account, and is related to the length of the three-phase boundary line. With perfect wetting ( $\theta=0^\circ$ ), this specific force reaches a maximum value, which corresponds to the surface tension of the liquid; with poor wettability ( $\theta>90^\circ$ ), however, this wetting force becomes negative. Energetic interactions between the solid and liquid phases are thus made visible through the contact angle. When the solid body is immersed, the *triple line scans the solid surface*, yielding the advancing contact angle ( $\theta_A$ ), and when it is withdrawn, the receding contact angle ( $\theta_R$ ) is obtained - measured by force of wetting. The difference ' $\theta_A-\theta_R$ ' is called the contact angle hysteresis (CAH).

**Physically:** The Wilhelmy equation  $\cos\theta_{(A|R)} = ^aF / (p \cdot \gamma)$  describes the relationship between the contact angle  $\theta$ , the perimeter of the solid  $p$ , the wetting force  $^aF$ , and the surface tension of the liquid  $\gamma$ . Adhesion tension  $^aE$  is defined according  $^aE_{A,R} = m_{meniscus} \cdot g / p = ^aF / p = \gamma \cdot \cos\theta = \gamma_s - \gamma_{sl}$

**Practically:** In this measurement, after the position of contact of liquid and solid has been determined, the 'Wilhelmy-Platte (DIN)' test specimen is moved into and out of the Isooktan 99.5+% liquid by incremental movements attended by 'CLT'(Constant-Level-Technic). This allows the determination of *most accurate static contact angles*. Furthermore, a defined force stabilization can be maintained before each subsequent movement. This enables *static contact angle* measurement based on the patent DE 4412405, which provides the often required equilibrium measurements in accordance with physical principles. This control technology yields a stabilization time for levelling and for each of the selected 94 measuring positions, the respective contact angle  $\theta_{A,s}$  or  $\theta_{R,s}$  is calculated independently.

**Information:** The IMETER M4 method description, available at '<https://imeter.de/?view=article&id=46>', provides information on procedures and calculations.

### 2. ===== Static Contact Angle =====

Wilhelmy-Platte (DIN), cubic plate 20×0,2mm


4,0mm Wilhelmy-Platte (DIN) \ Isooktan 99.5+%,  $\gamma = 18,75$  mN/m

$\theta = 20,0^\circ\text{C}$

Contact Angle, CA $\theta_c$	<0 ± 0°>		
Contact Angle Hysteresis, CAH <sub>s</sub>	<0°>		
	<b>Advancing <math>\theta_{A,s}</math></b>	<b>Receding <math>\theta_{R,s}</math></b>	<b>Transition - mm</b>
Contact Angles, $\theta_s$	0°	0°	0° <sub>adv.</sub> ⇒ 0° <sub>rec.</sub>
Rated measurements $n$ mm-range	41 0,000 - 4,000mm	53 4,000 - 0,000mm	2 4,000 - 3,975mm
avg. Triple line speed $\bar{v}_z$ [mm/s]	0,0286 ± 0,0115	-0,0241 ± 0,210	
avg. Step distance $\Delta h_0$ [mm]	0,100 ± 0	-0,0719 ± 0,0363	
avg. equilibration Time $\bar{t}_{eq}$ [s]	2,94 ± 0,386	2,86 ± 0,203	

*Static Contact Angle Measurement: For superwetting - consider adhesion parameters for analysis. The avg. equilibration Time  $t_{eq}$  is the mean of the time between positioning and acquiring the wetting force. The criterion of force equilibrium is defined in the utilized IMPro. Diagramm(s) below presenting these results and further details in a graphical way.*

→ Diagramm 1.2: 'Contact Angle<sup>2</sup>' -- Data, Calculation and Results

*Diagram 'Contact Angle<sup>2</sup>' gives a summary on this contact angle measurement: The small triangular marks show the individual contact angle values depending on the immersion depth of the sample. At '0' on the x-axis, the flatface front of the sample touches the surface of the liquid, to the right the depth of immersion increases. So 'immersion' indicates the nominal position of the three-phase boundary ('Triple Line') on the sample surface. Red triangles show the contact angles when the fluid advances over the sample surface, Green triangles when retracting (receding). The contact angle peak at 0,000 mm refers a sharp edge - the 'triple line' is fixed there while the immersion goes on until this hypertensive  $\theta_A$  max.  $0,0^\circ$  is reached. - The direction reversal shows normally (why not here?) a step, that is the 'contact angle hysteresis' (CAH). The two contact angles and the extent of hysteresis are determined from the graph by plotting linear regressions over the ranges. The advancing angle  $\theta_A$  is determined by the regression equation for the immersion depth '0'; the receding angle  $\theta_R$  is determined by extrapolation to the maximum immersion depth (thus any errors due to adherence of remaining liquid are eliminated). The standard deviation of the regression curves is indicated by the width of the hatching, which also graphically indicates*

the data range of the regression. Further marks indicate the determined angles (red, green) and the hysteresis range (yellow) above the Y-axis, the main result, the equilibrium contact angle  $\theta_c$  is shown in blue.

### Adhesion Force of the equilibrated Triple Line

Energy of Adhesion,  ${}^a\bar{E}_s = ({}^aE_A + {}^aE_R) / 2$  **18,887 ± 0,021 mN/m**

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
### 3. ===== Details on the Measurement and Setup =====

**Sample:** 'Wilhelmy-Platte (DIN)', **Form:** cubic plate, face dimension 20 x 0,2 mm.

**Test Liquid:** 'Isooktan 99.5+%' at 19,99°C: density 0,691916 g/mL, surface tension 18,75 mN/m, viscosity 0,496 mPa·s; Capillary length 1,7 mm.

**Suspension/Fixation:** by operator; measured weight at start 1,2118 g.

**Vessel:** tempered double wall measuring cell, height 140 mm, inclosing sample and liquid at coherent temperature and atmosphere;

Surface area  $\infty$  by CLT. - The 'CLT' Constant-Level-Technic prevents rising/falling of the Isooktan 99.5+% -level in the vessel (surface 1452mm<sup>2</sup>) through immersion/emersion of the sample volume in the vessel by appropriate pumping of Isooktan 99.5+%.

**Conditioning:** during 2,8 min in the measuring cell above the surface of the liquid.


**Sample Weight:** at start 1,5312 g, after conditioning 1,5311 g, after the measurement 1,6092 g.

*Loss in weight* during conditioning -0,0001 g (-33 ppm) but different after measurement 0,0780 g (+5,1%).

**Measurement algorithm:** static CA-measurement, acquisition of equilibrated values, stepwise movement. Maximum immersion of 4,000 mm and force equilibration at the inflection point for 3,3 s. Duration for immersion 2,3 min, for emersion 2,0 min.

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
3.3: 'Temperature & Events<sup>2</sup>' -- Survey on Temperature and Time

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- Diagramm 4.3: 'RawData<sup>21</sup>' -- Acquired weights during immersion and emersion (raw data and immersion corrected lengths)

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#### 4. ===== Table of raw data and results =====

The table below provides the data for each contact angle measuring position in this experiment ( $n=93$ ). - Within a row the column  $t$  lists the time of the CA-measurement.  $t_{eq.}$  is the duration reach the equilibrated force after positioning (static measurement). The nominal distance between the sample bottom flatface and the liquid surface is given by  $h_0$ , whilst  $h_{corr.}$  shows the corrected immersion depth- *due to CLT the values are the same*. With  $W_{RAW}$  the final total weight of the sample, sample holder, wetting force and buoyancy is written. The  $W_{corr.}$ -values\* are former weighing values when the weight of the sample and its holder is subtracted.

After buoyancy correction,  $\Delta F$  is the force exerted on the triple line. The Energy of Adhesion  $\Delta E_{A,R} = \Delta F/p$  is the force of adhesion per meter of the circumference ( $p$ ) of the sample.  $\theta$  is the Contact Angle. The last column indicates by **chr** the advancing resp. receding state.

(Resolving static force of wetting from  $\Delta F = (W_{RAW} - W_0 - W_{Buoy}) \cdot g + F_{Buoy,air}$ ) \* $W_{corr.} = W_{RAW} - W_0$  (with  $W_0 = 2,7429$  g).

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## Programm

Data created during execution of the IMPro "ContactAngle\_atConstantLevel", type 4/4. IMPro finished as projected.

(the digital twin of process IDN°20326 provides on demand reports with more details, audit-log protocol and structured raw data)

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